


# Disposable Feeding Needles



1. Special processed tips protect delicate organs
2. The straight tube allows smooth insertion into esophagus
3. EOG sterilized
4. Five needles per pouch, 20 pouches packed in a box

## For Mouse



Model	Size (mm)	Gauge	Material
493-18GS	φ 1.3 X L65	18G	Polyurethane



Model	Size (mm)	Gauge	Material
493-22CBS	φ 0.9 X L65	22G	Polyurethane



Model	Size (mm)	Gauge	Material
493-2230CBS	φ 0.9 X L30	22G	Polyurethane

## For Rat



Model	Size (mm)	Gauge	Material
493-16GLS	φ 1.7 X L80	16G	Polyurethane

V2-2404



**NATSUME SEISAKUSHO CO.,LTD.**

2-18-6, Yushima, Bunkyo-ku, Tokyo JAPAN 113-8551

<https://www.nazme.co.jp> / contact : [g-support@nazme.co.jp](mailto:g-support@nazme.co.jp)

